

Raman Spectroscopy for Detection of Critical Minerals and Rare Earth Elements (REEs) : Copper (Cu)

INTRODUCTION

Copper (Cu) is a critical industrial and strategic metal essential for electrical infrastructure, renewable energy systems, electric vehicles, power transmission, and electronics. With increasing global electrification and decarbonization efforts, demand for copper continues to rise, necessitating efficient mineral exploration and accurate identification methods.

Copper occurs in nature in various mineral forms, including carbonates, silicates, sulfides, and oxides. Reliable identification of copper-bearing minerals is vital for geological mapping, resource evaluation, and beneficiation planning.

Raman Spectroscopy has emerged as an effective analytical technique for copper mineral identification due to its ability to probe molecular and lattice vibrations that are highly specific to mineral structure.


MATERIALS AND METHODS

Raman spectra were acquired using the **TechnoS IndiRAM™ CTR Raman spectrometer**, designed to provide high spectral resolution, wavelength stability, and excellent signal-to-noise performance for mineralogical analysis.

Two structurally distinct copper-bearing minerals were Analyses:

- **Diopase** – $\text{Cu}_6\text{Si}_6\text{O}_{18} \cdot 6\text{H}_2\text{O}$ (hydrated copper cyclosilicate).
- **Malachite** – $\text{Cu}_2\text{CO}_3(\text{OH})_2$ (basic copper carbonate)

These minerals were selected due to their contrasting silicate and carbonate frameworks.



RESULTS AND DISCUSSION

Diopase

Diopase exhibits Raman features dominated by silicate ring vibrations:

- Cu–O lattice modes at ~ 137 , 165 , and 181 cm^{-1}
- SiO_4 bending vibrations in the $200\text{--}500\text{ cm}^{-1}$ range
- Internal silicate modes near 527 and 665 cm^{-1}
- Si–O symmetric stretching modes around 748 , 921 , 974 , and 1012 cm^{-1}

These sharp and well-defined peaks reflect the ordered cyclosilicate structure of Diopase.

Malachite

Malachite displays Raman signatures characteristic of carbonate minerals:

- Lattice vibrations at ~ 180 , 219 , 270 , 355 , 434 , 538 , and 596 cm^{-1} (with an intense Cu–O bending mode near $\sim 434\text{ cm}^{-1}$)
- Carbonate group vibrations, including:
- Symmetric bending at ~ 721 and 755 cm^{-1}
- Asymmetric bending near $817\text{--}820\text{ cm}^{-1}$
- Symmetric stretching near 1059 and 1097 cm^{-1}
- Asymmetric stretching modes around $1368\text{--}1490\text{ cm}^{-1}$

The dominance of carbonate-related bands clearly differentiates Malachite from silicate-based copper minerals.

Mineral Discrimination

Despite visual similarity, Diopase and Malachite show distinct and non-overlapping Raman fingerprints, enabling rapid and confident mineral identification without sample preparation.

CONCLUSION

Raman spectroscopy provides a fast, non-destructive, and highly specific method for identifying copper-bearing minerals. The distinct spectral signatures of Diopase (silicate) and Malachite (carbonate) demonstrate Raman's effectiveness for copper mineral discrimination in exploration and processing environments.

Using high-resolution Raman systems and with ongoing development of **Portable Raman solutions**, TechnoS Instruments enables accurate, field-ready mineral identification to support the growing demand for critical minerals such as copper.


REFERENCE

1. Rinaudo, C., Gastaldi, D., & Croce, G. Raman characterization of lithium aluminosilicates. *Journal of Raman Spectroscopy*, 36, 810–816 (2005).
2. Frost, R.L., et al. Raman spectroscopy of carbonate minerals. *Spectrochimica Acta Part A*, 67, 604–611 (2007).
3. McMillan, P.F. Vibrational spectroscopy of silicates. *Physics and Chemistry of Minerals*, 16, 245–254 (1988).
4. Data taken at IIGJ-Jaipur using IndiRAM CTR-300C Raman Spectrometer system.

OUR PRODUCTS



IndiRAM™ CTR Series



IndiRAM™ CTR for
Quantum Characterizer



IndiRAM™ CTR-
Mini Series



IndiRAM™ Edu




Diamond Detection
System




IndiRAM™
Handy Series



IndiRAM™ CTR
P-Series



IndiRAM™ Portable
System



G1-43, RIICO Industrial Area, Sitapura, Jaipur, Rajasthan
302022



www.technosphotronics.com



info@technos.in



+918003220052